Experiment 3: Reactions in Aqueous Solutions: Pre–lab

Due at the beginning of lab.

Name: \_\_\_\_\_\_\_
Section: \_\_\_\_\_\_/10

## 1. Precipitation Reactions

a. On the reverse side of this page or on a separate piece of paper, **neatly** write the <u>balanced</u> molecular equation, an ionic equation with spectator ions crossed out, and the balanced net ionic equation for the reaction of each pair of aqueous solutions. (Be sure to include all states, aq, l, s or g. Use the solubility rules as a guide.) Attach any extra pages with your work.

**Example:** lead nitrate and potassium iodide

Molecular equation:  $Pb(NO_3)_{2 \text{ (aq)}} + 2 \text{ KI}_{(aq)} \rightarrow PbI_{2 \text{ (s)}} + 2 \text{ KNO}_{3 \text{ (aq)}}$ 

Ionic equation:  $Pb^{2+}_{(aq)} + 2 NO_{3-}_{(aq)} + 2 K_{(aq)}^{+} + 2 I_{(aq)}^{-} \rightarrow PbI_{2(s)} + 2 K_{(aq)}^{+} + 2 NO_{3-}_{(aq)}$ 

Net ionic equation:  $Pb^{2+}_{(aq)} + 2 \Gamma_{(aq)} \rightarrow PbI_{2 (s)}$ 

- a) Sodium chloride and ammonium nitrate
- b) Sodium chloride and silver nitrate
- c) Sodium chloride and barium nitrate
- d) Sodium sulfate and ammonium nitrate
- e) Sodium sulfate and silver nitrate
- f) Sodium sulfate and barium nitrate
- g) Sodium hydroxide and ammonium nitrate

- h) Sodium hydroxide and silver nitrate
- i) Sodium hydroxide and barium nitrate
- j) Sodium carbonate and ammonium nitrate
- k) Sodium carbonate and silver nitrate
- I) Sodium carbonate and barium nitrate

b. Tabulate your expected results in the chart below based on the solubility rules covered in you text and lecture. Place a "P" in the box if you expect a precipitate to form and "NR" if you do not think a reaction will occur. Consider "slightly soluble" compounds as no reaction since they are closer to strong electrolytes than precipitates.

#### **Predicted Results**

	NH <sub>4</sub> NO <sub>3</sub>	AgNO <sub>3</sub>	Ba(NO <sub>3</sub> ) <sub>2</sub>
NaCl			
Na <sub>2</sub> SO <sub>4</sub>			
NaOH			
Na <sub>2</sub> CO <sub>3</sub>			

# **EXPERIMENT 3: REACTIONS IN AQUEOUS SOLUTIONS** (Read through this prior to beginning the pre—lab) This experiment corresponds to chapter 3 in your text)

#### Introduction

One of the most important characteristics of water is its ability to dissolve many compounds. Solutions in which water is the solvent are called *aqueous solutions*. Many important reactions take place in aqueous solutions. In fact, many of the reactions that take place throughout your body (from your organs down to individual cells) are aqueous reactions. Understanding the most common aqueous reactions and how to correctly write them is one of the most important skills you should master in Chemistry 1A. This skill will be used extensively throughout the remainder of the semester and in Chemistry 1B.

Before you begin, one must ask the question: What observations indicate that a chemical reaction has occurred? Some indications include: the formation of an insoluble solid (precipitate), color change, the evolution of a gas, or a temperature change. In this experiment, you will predict what will happen when two aqueous solutions are mixed, and then test your predictions in the laboratory. During the previous discussion period, your lab instructor lectured on the topic of reactions in aqueous solution with examples of the correct way to write a molecular equation, an ionic equation, and the overall net ionic equation for several types of aqueous reactions.

In exchange or double displacement reactions, cations and anions exchange partners as in the following generic reaction:  $AX + BY \rightarrow AY + BX$ 

There are three types of exchange reactions—precipitation reactions, gas-forming reactions, and neutralization reactions. An example of each type is given below:

Reaction Type	<u>Example</u>
Precipitation—results in the formation of an insoluble solid called a <i>precipitate</i>	$Pb(NO_3)_2$ (aq) + 2 KI (aq) $\rightarrow$ 2 KNO <sub>3</sub> (aq) + $PbI_2$ (s) precipitate
Gas-forming—results in the formation of a gas that escapes from the solution	2 HCl (aq) + Na <sub>2</sub> S (aq) $\rightarrow$ H <sub>2</sub> S (g) + 2 NaCl (aq)
Neutralization-reaction of an acid and a base which results in the formation of a salt and water	$HNO_3$ (aq) + KOH (aq) $\rightarrow$ KNO <sub>3</sub> (aq) + H <sub>2</sub> O (/) acid base salt water

A exchange reaction will occur if (1) a precipitate forms from soluble reactants or (2) a stable molecule forms, such as water or an insoluble gas. Refer to the example on the next page and the sections on ionic equations and exchange reactions in your textbook for information on writing chemical equations for exchange reactions. Familiarity with the solubility rules tabulated on the next page is required to write these equations.

**Objectives:** Upon completion of this exercise and laboratory experiment, you should be able to:

- 1. Using the solubility rules, determine the species present in aqueous solutions of compounds.
- 2. Predict the type of reaction that will occur when two aqueous solutions are mixed.
- 3. Write the chemical equation, the ionic equation, and the net ionic equation for reactions taking place between aqueous solutions.
- 4. Experimentally identify the type of reaction occurring when two aqueous solutions are mixed through simple visual and temperature measurements.

## **EXPERIMENT 3: REACTIONS IN AQUEOUS SOLUTIONS:**

Solubility Rules: You should be familiar with these rules and be able to use them efficiently.

Soluble Compounds Compounds containing NO <sub>3</sub> <sup>-</sup> and C <sub>2</sub> H <sub>3</sub> O <sub>2</sub> <sup>-</sup> ClO <sub>3</sub> <sup>-</sup> , ClO <sub>4</sub> <sup>-</sup> , Gr I cations and NH <sub>4</sub> <sup>+</sup>	Exceptions None
Compounds containing Cl <sup>-</sup> , Br <sup>-</sup> , and l <sup>-</sup>	Compounds containing Pb <sup>2+</sup> , Ag <sup>+</sup> & Hg <sub>2</sub> <sup>2+</sup>
Compounds containing $SO_4^{\ 2-}$	Compounds containing $Ca^{2+}$ , $Hg_2^{2+}\{ss\}$ , $Ag^{+}\{ss\}$ , $^1$ $Sr^{2+}$ , $Ba^{2+}$ , and $Pb^{2+}$
Insoluble Compounds	<b>Exceptions</b>
Compounds containing CO <sub>3</sub> <sup>2-</sup>	Compounds containing NH <sub>4</sub> <sup>+</sup> and Gr I metal cations
Compounds containing OH <sup>-</sup>	Compounds containing NH <sub>4</sub> <sup>+</sup> , Ca <sup>2+</sup> {ss}, Sr <sup>2+</sup> {ss}, Ba <sup>2+</sup> {ss}, and Gr I metal cations
Compounds containing S <sup>2-</sup>	Compounds containing $\mathrm{NH_4}^+$ and $\mathrm{Gr}\ \mathrm{I}\ \mathrm{metal}\ \mathrm{cations}$

In your reading you may encounter salts that are classified as "slightly soluble". For the purposes of this course, we will consider them to be completely soluble. Later in chem. 1B, you will understand the extent of solubility quantitatively when you learn more about equilibrium.

Writing the molecular, ionic, and net ionic equations for a chemical reaction in aqueous solution:

The steps for writing a net ionic equation for the exchange reaction of aqueous  $K_2SO_4$  and  $Pb(C_2H_3O_2)_2$  are outlined below.

1. Write the balanced *molecular equation* which gives the complete chemical formula and the phase of the reactants and products. Switch the reactant cation/anion pair to form the products and use the solubility rules to determine the phase of the products.

$$K_2SO_{4 (aq)} + Pb(C_2H_3O_2)_{2 (aq)} \rightarrow PbSO_{4 (s)} + 2 KC_2H_3O_{2 (aq)}$$

2. Write the balanced *complete ionic equation* which shows all strong electrolytes as ions.

$$2 K^{+}_{(aq)} + SO_{4}^{2-}_{(aq)} + Pb^{2+}_{(aq)} + 2 C_{2}H_{3}O_{2}^{-}_{(aq)} \rightarrow PbSO_{4 (s)} + 2 K^{-}_{(aq)} + 2 C_{2}H_{3}O_{2}^{-}_{(aq)}$$

3. Cancel out the *spectator ions*, the ions that appear as both reactants and products in the complete ionic equation. What remains is the *net ionic equation*.

$$\frac{2 K^{+}_{(aq)} + SO_{4}^{2-}_{(aq)} + Pb^{2+}_{(aq)} + 2 C_{2}H_{3}O_{2}^{-}_{(aq)}}{Pb^{2+}_{(aq)} + SO_{4}^{2-}_{(aq)}} \rightarrow PbSO_{4}_{(s)} + \frac{2 K^{+}_{(aq)} + 2 C_{2}H_{3}O_{2}^{-}_{(aq)}}{Pb^{2+}_{(aq)} + SO_{4}^{2-}_{(aq)}} \rightarrow PbSO_{4}_{(s)}$$

<sup>1</sup> Compounds that are "slightly soluble" {ss}, are considered to be soluble in this course. These compounds are closer to strong electrolytes than precipitates in terms of solubility.

Adapted from "Reactions in Aqueous Solutions" by David Reichgott and Mary O'Brien, Edmonds Community College, Lynnwood, Washington and "Reactions in Aqueous Solutions" Illinois State University, Normal, Illinois.

# Experiment 3: Reactions in Aqueous Solutions: Experimental Procedure:

In this part of the experiment, you will compare you pre—lab assignment predictions of reactions in aqueous solutions against your observed results obtained when aqueous solutions are mixed in the laboratory.

- 1. Unless noted, the test tubes that you use do not need to be completely dry. However, rinse them with deionized water before use.
- 2. Select three small test tubes. Using a graduated cylinder, measure 1.0 mL of deionized water and transfer it to a test tube. Use a grease pencil to mark the water level in the test tube and then use the test tube as a guide to make a 1.0 mL mark on the other two test tubes. Repeat this step to make 2.0 mL marks on each test tube. Each of the three test tubes should have a 1.0 mL and a 2.0 mL mark when you are finished.

### **Precipitation Reactions:**

1. Add 1 mL of 0.1 M NaCl to each of the three small test tubes. Then, in order, add 1 mL of each of the solutions listed in the first row of the table in question 1 of your lab report. For example, the first test tube should contain NaCl and NH<sub>4</sub>NO<sub>3</sub>, the second should contain NaCl and AgNO<sub>3</sub>, and the third should contain NaCl and Ba(NO<sub>3</sub>)<sub>2</sub>. Record your observations. Dispose of the solutions in the appropriate waste container.

Repeat step 1 for all solution combinations in the table. Again, dispose of all solutions in the appropriate waste container.

### **Gas-Forming Reactions:**

- 2. To a <u>dry</u> test tube, add a sample of CaCO<sub>3</sub> that is about the size of a pencil eraser. Then add 1 mL of deionized water followed by 1 mL of 2 M HCl. Record your observations.
- 3. Repeat step 2 using Na<sub>2</sub>SO<sub>4</sub> in place of CaCO<sub>3</sub>.

#### **Neutralization Reactions:**

- 4. Add a scoop of about 0.5 g of Mg(OH)<sub>2</sub> (about the size of your pinky nail) into a 100 ml beaker. Use a graduated cylinder to measure 25.0 mL of deionized water. Pour the water into the beaker containing Mg(OH)<sub>2</sub>. Swirl the mixture.
- Use a graduated cylinder to measure 25.0 mL of 1.0 M HCl. Add the HCl to the beaker containing the Mg(OH)<sub>2</sub> mixture and thoroughly stir the solution. Observe what happens, upon mixing.

EXP 3:	REACTIONS IN AC	QUEOUS SOLUTIONS	Name:		
EXPERIMENTAL RESULTS		Section:	Score:	/ 30	
-	ions. Use an asteri	Tabulate your observationsk (*) to mark any obser			
		NH <sub>4</sub> NO <sub>3</sub>	AgNO <sub>3</sub>	Ba(NO <sub>3</sub> ) <sub>2</sub>	
NaCl					
Na <sub>2</sub> SO	1				
NaOH					
Na <sub>2</sub> CO	3				
2. 003-	molecular equation	ppened when you mixed on for the reaction.  ppened when you mixed on for the reaction.			
3.		eactions: when you mixed the Mg( e molecular equation fo		P Explain your observat	ions. If there was a
		Instructo	r Date and Sign:		

Experiment 3

Chem.1A

CSUS Department of Chemistry

EXP 3:	REACTIONS IN AQUEOUS SOLUTION	ONS	Name:			
Predict	ions		Section:			
1.	Precipitation Reactions: Comple	eted in p	art 1 of this assig	nment.		
2.	Gas-Forming Reactions (refer to )	Gas-Forming Reactions (refer to you text and notes)				
	Which of the following pairs of re	eactants	will give off a gas	when mixed?		
	V	Will a gas	s evolve? (Y/N)	Formula of (	Gas	
	HCl (aq) and CaCO <sub>3</sub> (s)				_	
	HCl (aq) and Na₂SO₄ (aq)				_	
ionic ed	Neutralization Reactions he balanced molecular equation, a quation for the reactions below. (B ations of sodium hydroxide and hyd	Be sure to	o include all state		d out, and the bala	nced net
(b) Solu	itions of hydrofluoric acid and calc	cium hyd	lroxide.			
(c) Solu	tions of nitric acid and aqueous ar	mmonia.				
(d) Soli	d magnesium hydroxide and aqued	ous perc	chloric acid.			

Experiment 3

Chem.1A

CSUS Department of Chemistry

Classification:

c. Cadmium chloride and sodium sulfide

Classification: