

Result of Standard chromatogram by Instructor			
Compound	Retention time (min)	Peak Area ( $\mu\text{s} \times \text{min}$ )	Conc. (ppm)
Fluoride	2.75	0.976	4.88
Chloride	3.70	3.836	20.0
Nitrate	5.72	3.858	40.0
Sulfate	6.72	5.207	40.0

Chromatogram of my Solution A (IC # 1)		
Compound	Retention time (min)	Peak Area ( $\mu\text{s} \times \text{min}$ )
Unknown	5.80	3.427
Compare to Standard compounds retention time: my Unknown is : <b>NO<sub>3</sub><sup>-</sup></b>		

Estimated concentration of unknown:

$$\frac{A_{unk}}{A_{known}} = \frac{C_{unk}}{C_{known}}$$

$$\frac{3.427}{3.858} = \frac{C_{unk}}{40.0}$$

$$C_{unk} = 35.5 \text{ ppm}$$

Note: Do same calculation for the concentration of  $[\text{F}^-]$  in your home's tap water.

Make two lower and higher standard solutions (use volumetric flask & pipet):

[50% of conc. Of unknown]:

$$35.5 \text{ ppm} \times 0.5 = 17.75 \text{ ppm}$$

Stock solution  $[\text{NaNO}_3] = 2000 \text{ ppm}$

$$2000 \times V = 17.75 \times 100$$

$$V = 0.8875 \text{ mL} \quad (\text{use } 1.0 \text{ mL volumetric pipet})$$

Re-calculate for the new conc.:

$$2000 \times 1.0 = C \times 100$$

$$C = 20 \text{ ppm}$$

[150% of conc. Of unknown]:

$$35.5 \text{ ppm} \times 1.50 = 53.25 \text{ ppm}$$

Stock solution  $[\text{NaNO}_3] = 2000 \text{ ppm}$

$$2000 \times V = 53.25 \times 100$$

$$V = 2.6625 \text{ mL} \quad (\text{use } 3.0 \text{ mL volumetric pipet})$$


Re-calculate for the new conc.:

$$2000 \times 3.0 = C \times 100$$

$$C = 60 \text{ ppm}$$

Quantitative Analysis: Running Solutions A , STD 1 and STD 2		
Solution	Conc. (ppm)	Peak Area ( $\mu\text{s} \times \text{min}$ )
A	?	3.580
STD1	20	1.824
STD2	60	5.539
Instructor STD	40	3.858

Plot peak area (A, STD1, STD2 and Instructor STD) versus concentration of nitrate ion:



Use the line equation from your graph:

$$P_{\text{area}} = 0.0929 C + 0.0253$$

Use the peak area of the solution A from the quantitative chromatogram:

$$3.580 = 0.0929 C + 0.0253$$

$$C = 38.2 \text{ ppm}$$