

Lab 7: pAMP & pKAN ligation

Objective of this lab is to create a **recombinant plasmid** constructed from pieces of two different plasmids joined together.

Define: recombinant / cloning / vector / insert

This will be done by:

- Digesting both plasmids with the same restriction enzymes to create compatible **sticky ends**
- Completeness of enzymatic digestion will be evaluated by running a small aliquot on agarose

This is crucial because uncut plasmids transform at a MUCH HIGHER efficiency than recombinant ones.

If uncut plasmid is included in the ligation & transformation, the majority of transformants will contain the *original, uncut plasmids* NOT the recombinant.

Double transformants (which contain both original, uncut plasmids) will be phenotypically identical to desired clone

- Restriction enzymes will be **heat inactivated** to prevent competition with the ligation
- Cut plasmids will be mixed and sticky ends allowed to **anneal**
- **DNA ligase** will be used to covalently link the annealed fragments together
 - DNA ligase* uses the energy of an ATP molecule (and Mg^{++} cofactor) to restore the phosphodiester bond between the 5' phosphate and 3' OH on the DNA fragments (in other words, **it seals nicks**).
 - Sticky end ligations are performed at room temperature, a compromise between enzymatic efficiency (best at 37C) and annealing of sticky ends (best in the cold)
- Because we are not specifically isolating the two fragments we wish to ligate together, a variety of products of ligation are possible. The desired product (called pAK) must have the following:
 1. An **origin of replication (ori)**
 2. The **amp** gene for ampicillin resistance
 3. The **kan** gene for kanamycin resistance
- Transformants containing this recombinant plasmid will be dual resistant (able to grow on agar containing both ampicillin and kanamycin)

Ampicillin inhibits cell wall synthesis, **bacteriostatic**

Kanamycin, an aminoglycoside antibiotic, inhibits protein expression (30S subunit of the 70S ribosome), **bactericidal**

- As a consequence, the *recovery step* following heat shock during transformation is crucial


Blunt end ligations are also possible (though we are not doing any here).

- DNA is cut by restriction enzymes that do *not* leave a single strand overhang (sticky end)
- Or, sticky ends are digested away by **single-strand specific nuclease** such as mung bean nuclease, leaving blunt ends
- Blunt end ligations are much less efficient than sticky ones:
 - contact between DNA fragments is transient; no hydrogen bonds hold them together.
 - Nonspecific joining of any fragment to any other fragment

Cloning strategies & issues


Example #1.

You want to clone the human insulin gene into a plasmid vector so you can sequence it.



Cut the insulin gene out of chromosome 11 using HindIII. Purify this fragment. Open the plasmid pGEM at its multicloning site using HindIII. Mix and ligate; transform bacteria; select for bacteria that can grow on ampicillin.

You get TONS of colonies/transformants. You miniprep 4 of them, and cut the plasmid DNA with HindIII. Run a gel.



What happened? What would a digest of the desired plasmid look like?

Phosphatase enzymes remove exposed phosphate groups from the ends of DNA.

Remember: DNA ligase requires that the DNA pieces it is joining together have a 5' phosphate to form the new phosphodiester bond.


Solution: As before, cut the pGEM vector with HindIII. Next, treat it with a phosphatase (CIP, or calf intestinal phosphatase, is commonly used). This will remove the 5' phosphates and make it impossible for the vector to re-close on itself during the ligation.

Example #2.

Same as above, except insulin gene is now flanked by EcoRI sites, not HindIII.

How can you clone an EcoRI fragment into a HindIII site?

Solution: By **blunt end ligation**.



To change a sticky end to a blunt end, you must remove the single strand overhang. This can be done by a **single strand specific nuclease** (a popular choice is mung bean nuclease). It chews up the single stranded part of the DNA end, and stops when it hits the double stranded part.


IMPORTANT: Creating a blunt end this way **destroys the enzyme recognition site**. After cloning, you cannot pop your insert out again with EcoRI *or* HindIII.

Blunt end ligation is possible. DNA ligase can link blunt ends together just as it can link sticky ends; it is less efficient because the blunt ends are not held together by any hydrogen bonds.

To compensate for this inefficiency, **blunt end ligations** will include a **higher insert:vector ratio** in the ligation mix. By having, say, 10 molecules of your insert DNA for every one molecule of vector, you greatly increase the likelihood that any plasmid being transformed from this ligation will actually include an insert.

Example #3. Directional cloning


You want to *express* the human insulin gene inside a eukaryotic cell.



This is an example of a plasmid vector that will direct expression of a cloned gene in eukaryotic cells.

KEY POINTS:

- **Promoter** drives expression of genes downstream in a specific direction (see arrow)
- **Orientation** of inserted gene in multicloning site is crucial (gene is “sense” in only one direction)




How can you make this clone?

Cloning techniques & strategies

Why use two different restriction enzymes?

Using 2 enzymes to generate fragments with incompatible sticky/cohesive ends greatly improves your chances of obtaining the desired product of ligation (compared to using a single restriction enzyme). Consider the drawing below:




You want to move a cloned gene of interest from pGENE into pEMPTY. You need GENE inserted into the EcoRI site of pEMPTY to place it under the control of a promoter upstream from that site. You cut pGENE and pEMPTY with EcoRI, mix, and ligate.

◆ Site of action of ligase; requires energy

Predicting Products of a Ligation

Shown below are the fragments generated by complete digestion of pAMP and pKAN with BamHI and HindIII. Draw the listed ligation products and identify which ones, when transformed into competent bacteria, will “grow” on amp⁺ plates, kan⁺ plates, and amp⁺/kan⁺ plates.



Recognize that multimers of more than 2 ligated fragments are rare unless enzymatic digestion was incomplete. This is because it is a rare event for two DNA fragments, floating freely in solution, to find each other at all; it is even more rare for several to bind together before being ligated shut.

Possible products of ligation (indicate which is the desired product):

1+2	1+3	1+4	2+3
2+4	3+4		

Growth?
Amp+

Kan+

Amp/Kan+

1+1	2+2	3+3	4+4
1+2+3+4			

Growth?
Amp+

Kan+

Amp/Kan+