

Instructions for importing data into Excel for experiment 14

1. Open Excel.
2. Click on the **Data** tab.
3. Click Get **External Data** under the Data Tab.
4. Click on **Import Text File**. Select file on USB key.
5. Text Import Wizard pops up: You should see data from your file in the window.
6. Make sure **Delimited** is selected, Not Fixed width.
7. Click **Next**.
8. Click **Comma** box and leave Tab box checked. Click **Next**.
9. Click **Finish**.
10. Your data should appear in the spreadsheet with wavelength (nm) on the left followed by absorbance on the right.
11. Rearrange data as needed and plot your graph.